(en) ALYX 2RBC LR Kit**For Collection of Two Leukoreduced Red Blood Cell Products by Centrifugation.
For Use With the ALYX System.**

Store at Controlled Room Temperature. Protect from freezing. Avoid excessive heat.

Definition of "Controlled Room Temperature":

"A temperature maintained thermostatically that encompasses the usual and customary working environment of 20° to 25°C (68° to 77°F); that results in a mean kinetic temperature calculated to be not more than 25°C; and that allows for excursions between 15°C and 30°C (59° and 86°F) that are experienced in pharmacies, hospitals, and warehouses. Provided the mean kinetic temperature remains in the allowed range, transient spikes up to 40°C are permitted as long as they do not exceed 24 hours ... The mean kinetic temperature is a calculated value that may be used as an isothermal storage temperature that simulates the non isothermal effects of storage temperature variations."

Reference: United States Pharmacopeia, General Notices. United States Pharmacopeial Convention, Inc. 12601 Twinbrook Parkway, Rockville, MD.

Description:

The ALYX 2RBC-LR Kit is a single unit consisting of:

- (1) separation chamber
- (1) umbilicus
- (1) cassette with return line filter with pore size of 170 microns
- (1) 0.9% sodium chloride solution container
- (1) ACD-A anticoagulant solution container
- (1) ADSOL red cell preservation solution container
- (3) blood component containers
- (2) red cell storage containers
- (1) 18 gauge needle
- (1) pre-donation sampling pouch
- (1) Asahi Kasei Medical Co. SEPACELL leukoreduction filter

Indications and Usage:

The system is an automated apheresis system designed to collect red blood cells (RBCs) by centrifugation. The RBCs must be handled per center requirements and applicable local regulations.

If ACD-A/ADSOL red cells are to be frozen, they must be frozen within 6 days from date of collection.

ACD-A/ADSOL red cells may be irradiated. If irradiation dose is 3000 cGy, red cells should be transfused within 28 days from date of collection.

Contraindications:

Use of the ALYX component collection system is contraindicated when adequate anticoagulation cannot be achieved.

Warnings:

As with any procedure involving extracorporeal circulation, there is a risk of hemolysis, air embolism or blood clotting. To avoid these risks, do not deviate from the Operator's Manual for the ALYX component collection system.

Approximate Container Tare Weights:

Red Cell container with label	31 g
In-Process container	30 g
Plasma container	30 g
Red Cell storage container (1) with label, segment tubing, 6 inches of kit tubing, and no clamp.	43 g

Nominal Tubing Dimensions:

	ID	OD
Kit tubing - excluding needle tubing	0.118 inches 3.0 mm	0.178 inches 4.5 mm
Needle tubing	0.130 inches 3.3 mm	0.205 inches 5.2 mm
Segment tubing	0.126 inches 3.2 mm	0.178 inches 4.5 mm

Approximately 0.2 mL per inch (0.08 mL per cm).

Precautions:


1. Allow kit to reach room temperature before use.
2. Before operating the ALYX system, the operator should be trained and knowledgeable in the system components, controls and operating procedure.
3. If closed system conditions are maintained, the red cell product may be stored at 1 - 6°C for up to 42 days. If open system collection conditions occur, the red cell product must be transfused within 24 hours of collection.
4. Red blood cells collected from some donors may have extended filtration times resulting in the potential for ineffective filtration and leukoreduction. The final product (after filtration) should be evaluated for residual leukocytes and red cell quality any time a reduction of flow through the leukoreduction filter occurs during the automated leukoreduction filtration process. The procedure information screen will display a message if the product needs to be checked.
5. Do not strip the tubing below the leukoreduction filter until separated from the filter.
6. Failure to follow instructions prior to venipuncture may result in dilution of whole blood sample.
7. Position the sample pouch to prevent air from entering the blood collection tube(s).

Directions for Use:

Refer to the ALYX component collection system Operator's Manual for the installation procedure for the disposable kit and complete instructions for use.

Dispose of waste in appropriate biohazard container or according to local regulatory requirements.

Actual kit may vary slightly from this diagram.

Symbols with Definitions:

	Caution, consult instructions for use
	Sterilized by a combination of steam and radiation. Sterile fluid path.
	Non-pyrogenic fluid path
	Do not reuse
	Do not use if the product sterile barrier system is compromised
	Do not vent
	Code
	Lot
	Expiration date
	Manufacturer
	Manufacturing facility / Manufactured by

	Fragile
	This way up
	Recyclable
	Red Cell Preservation Solution Container
	Saline Container
	Anticoagulant Solution Container
	Plasma container
	Red cell container
	In-process container

Rx only

For US Only. United States federal law restricts this device to sale by or on the order of a licensed health care practitioner.

SEPACELL is a trademark of Asahi Kasei Medical Co., Ltd.

© 2022 Fresenius Kabi AG. All rights reserved.

Fresenius Kabi AG
61346 Bad Homburg / Germany
Tel.: +49 (0) 61 72 / 686-0
www.fresenius-kabi.com

Fenwal International, Inc.
Camino Real Industrial Park,
Road No. 122
San German, Puerto Rico 00683

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Referencia: United States Pharmacopeia, Avisos generales. United States Pharmacopeial Convention, Inc. 12601 Twinbrook Parkway, Rockville, MD.

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© 2022 Fresenius Kabi AG. Todos los derechos reservados.
Fresenius Kabi México S.A. de C.V.
Fresenius Kabi S.A.